

# **OPENFDA MEDICAL DEVICE ADVERSE EVENT ANALYSIS**

Analysis of FDA Medical Device Adverse Event  
Reports

# Table of Contents

<b>1. Executive Summary</b>	<b>2</b>
<b>2. Methodology</b>	<b>2</b>
2.1. Data Source . . . . .	2
2.2. Data Standardization . . . . .	2
2.3. Problem Classification . . . . .	2
<b>3. Temporal Trend Analysis</b>	<b>3</b>
3.1. Overall Reporting Trends . . . . .	3
3.2. Cumulative Reports Over Time . . . . .	3
<b>4. Product Problem Analysis</b>	<b>5</b>
4.1. Top Product Problems . . . . .	5
<b>5. Patient Problem Analysis</b>	<b>7</b>
5.1. Patient Problems Overview . . . . .	7
<b>6. Manufacturer Analysis</b>	<b>8</b>
6.1. Top Manufacturers . . . . .	8
<b>7. Device Brand Analysis</b>	<b>9</b>
7.1. Top Device Brands . . . . .	9
7.2. Brand Temporal Trends . . . . .	11
<b>8. Technical Appendix</b>	<b>11</b>
8.1. Data Source . . . . .	11
8.2. Analysis Tools . . . . .	12
8.3. Exclusion Criteria . . . . .	12
8.4. Report Metadata . . . . .	12

# 1. Executive Summary

This report analyzes **4,262** FDA medical device adverse event reports submitted between **January 01, 2025** and **June 01, 2025** (a period of **5** months). The dataset includes **16** unique manufacturers and **63** unique device brands. The average reporting rate was **710.3** reports per month, with peak reporting of **1208** reports in **April 2025**.

Table 1: Summary statistics of FDA MAUDE reports

Metric	Value
Total Reports	4,262
Date Range	January 01, 2025 to June 01, 2025
Reporting Duration	151 days ( 5 months)
Unique Manufacturers	16
Unique Device Brands	63
Average Monthly Reports	710.3 reports/month
Maximum Monthly Reports	1208 reports in April 2025

# 2. Methodology

## 2.1. Data Source

The search of the FDA Manufacturer and User Facility Device Experience (MAUDE) database was conducted using the openFDA Device Event API. Retrieved JSON-formatted data were downloaded directly from the API and subsequently processed and cleaned for analysis.

## 2.2. Data Standardization

The analysis employs fuzzy matching algorithms to standardize manufacturer and brand names, addressing inconsistencies in naming conventions across reports. This standardization process uses the **RapidFuzz** library with partial ratio matching to group similar names under a canonical representation.

## 2.3. Problem Classification

Patient and product problems are extracted and categorized, excluding non-informative categories to focus on meaningful data patterns.

## 3. Temporal Trend Analysis

### 3.1. Overall Reporting Trends




Figure 1: Monthly trend of FDA MAUDE reports

#### **i** Trend Observation

The reporting shows variability across months, with notable peaks in **April 2025, March 2025**. The average monthly reporting rate is **710.3** reports, with a standard deviation of **393.4** reports.

### 3.2. Cumulative Reports Over Time




Figure 2: Cumulative FDA MAUDE reports over time

# 4. Product Problem Analysis

## 4.1. Top Product Problems




Figure 3: Most frequently reported product problems

### Product Issues

The most frequently reported problems are:

1. **Break** - 407 reports (9.5%)
2. **Scratched Material** - 324 reports (7.6%)
3. **Material Discolored** - 286 reports (6.7%)

These top three issues account for **1017** reports (**23.9%** of all reports).

# 5. Patient Problem Analysis

## 5.1. Patient Problems Overview




Figure 4: Patient problems reported in FDA MAUDE events



## 6. Manufacturer Analysis

### 6.1. Top Manufacturers




Figure 5: Top 10 manufacturers by report volume

The top manufacturer, **ALCON LENSX INC**, accounts for **1981** reports (**46.5%** of total).

Table 2: Top 5 manufacturers with report statistics

Rank	Manufacturer	Reports	% of Total
1	ALCON LENSX INC	1,981	46.48%
2	ALCON LABORATORIES IRELAND LTD	1,144	26.84%
3	STAAR SURGICAL	710	16.66%
4	RAYNER INTRAOCULAR LENSES LIMITED	145	3.40%
5	BAUSCH LOMB	107	2.51%

## 7. Device Brand Analysis

### 7.1. Top Device Brands




Figure 6: Top 10 device brands by report volume

The top device brand, **CLAREON TORIC IOL**, accounts for **910** reports (**21.4%** of total).

Table 3: Top 5 device brands with report statistics

Rank	Brand	Reports	% of Total
1	CLAREON TORIC IOL	910	21.35%
2	IMPLANTABLE COLLAMER LENS	710	16.66%
3	AKREOS INTRAOCULAR LENS	356	8.35%
4	MONARCH IOL DELIVERY SYSTEM CARTRIDGE	315	7.39%
5	CLAREON IOL	308	7.23%

## 7.2. Brand Temporal Trends




Figure 7: Cumulative adverse event reports by top 5 brands

## 8. Technical Appendix

### 8.1. Data Source

- **Database:** OpenFDA Medical Device Adverse Events
- **API Access:** <https://open.fda.gov/apis/device/event/>

## 8.2. Analysis Tools

- **Python 3.x** with pandas and rapidfuzz for data processing
- **R** with tidyverse, ggplot2, and lubridate for visualization
- **Fuzzy Matching:** RapidFuzz library for text standardization
- **Visualization:** ggplot2 with custom base R theme
- **Report Generation:** Quarto for reproducible analytics

## 8.3. Exclusion Criteria

### Patient Problems Excluded:

No Code Available, No Known Impact, Symptoms or Conditions,  
No Information, No Consequences, Insufficient Information

### Product Problems Excluded:

Adverse Event Without Identified Device, No Apparent Adverse Event,  
Appropriate Term/Code Not Available, Unknown, Insufficient Information

## 8.4. Report Metadata

- **Generated:** 2025-11-19 14:41:20
- **Dataset Version:** 2025-06-01
- **Total Records Analyzed:** 4,262
- **Analysis Pipeline:** OpenFDA Medical Device Adverse Event Analysis v2.1

## List of Figures

Figure 1	Monthly trend of FDA MAUDE reports .....	3
Figure 2	Cumulative FDA MAUDE reports over time .....	4
Figure 3	Most frequently reported product problems .....	5
Figure 4	Patient problems reported in FDA MAUDE events .....	7
Figure 5	Top 10 manufacturers by report volume .....	8
Figure 6	Top 10 device brands by report volume .....	10
Figure 7	Cumulative adverse event reports by top 5 brands .....	11

## List of Tables

Table 1	Summary statistics of FDA MAUDE reports .....	2
Table 2	Top 5 manufacturers with report statistics .....	8
Table 3	Top 5 device brands with report statistics .....	11